


# inada CUBE

FML-3000A



## CONTENTS

■ SAFETY PRECAUTIONS.....	1
■ COMPONENTS & DESCRIPTION .....	7
■ MAIN FEATURES .....	9
■ BEFORE USE .....	10
■ INTENSITY ADJUSTMENT .....	14
■ OPERATION .....	15
■ AFTER USE .....	16
■ TROUBLESHOOTING.....	17
■ FAQ .....	18
■ CARE AND STORAGE .....	19
■ CUSTOMER SERVICE AND WARRANTY .....	21
■ REGULAR MAINTENANCE FOR SAFETY .....	22
■ SPECIFICATIONS.....	<b>Back cover</b>

### Operating Manual



Warranty card provided separately

- Please read this Operating Manual before using this product.
- Keep this manual readily available for your reference.
- This product is designed for home use only.





# SAFETY PRECAUTIONS

Be sure that you have read these SAFETY PRECAUTIONS thoroughly before using your INADA CUBE for the first time. The Operating Manual contains important information about how to use the chair safely and correctly, avoiding possible damage to property or injury to persons using the chair and others. The Symbols shown below are used throughout the Operating Manual to highlight important safety instructions. Please adhere to all precautions. Descriptions of each component can be found on pages 7-8.

## Display Warnings




 <b>Warning</b>	This action could result in serious injury or death.
 <b>Caution</b>	This action could result in personal injury or damage to property.

## Types of Symbols and Description





	 Prohibited (The symbol on the left means "Do not disassemble.")
	 Required Action (The symbol on the left means "Unplug the chair.")

Remember to keep this Operating Manual in an easily accessible place.

## Safety Precautions: Installation

 <b>Caution</b>	
	<ul style="list-style-type: none"> <li>● <b>Do not install the chair in an area of high humidity, such as a bathroom or sauna.</b> → Could result in electric shock or damage to the chair.</li> </ul>
	<ul style="list-style-type: none"> <li>● <b>Only install on a flat, level surface.</b> → Chair could tip over accidentally if placed on an uneven surface.</li> </ul>

## Safety Precautions: Before Use

 <b>Warning</b>	
	<ul style="list-style-type: none"> <li>● <b>The following persons must not use this massage chair:</b> <ul style="list-style-type: none"> <li>○ Persons who have been instructed by a doctor to avoid massage therapy. Ex. Thrombosis, severe aneurism, acute varicosity, or dermatitis and/or skin infection including inflammation of the typodermis.</li> <li>○ Persons with suspected acute neck sprain (whiplash injury)</li> </ul> </li> </ul>
	<ul style="list-style-type: none"> <li>● <b>Persons with any of the following conditions should seek medical advice before using the chair.</b> <ul style="list-style-type: none"> <li>○ Serious heart conditions</li> <li>○ Sensory impairment caused by serious peripheral circulatory disorders due to diabetes and other diseases</li> <li>○ Osteoporosis, spinal fractures, acute pain and/or symptoms thereof</li> <li>○ Persons being treated for personal injuries</li> <li>○ High fever (100°F or higher)</li> <li>Ex. Strong acute inflammation illness caused by fatigue, chills, blood pressure fluctuation, or asthenia</li> <li>○ Persons in the early stage of pregnancy, or immediately after giving birth</li> <li>○ Persons who have electronic medical device(s) easily susceptible to electromagnetic interference</li> <li>Ex. pacemaker</li> <li>○ Abnormal or curved spine</li> <li>○ Persons planning to use the chair for treatment or for rehabilitation purposes</li> <li>○ Persons with malignant tumor(s)</li> <li>○ Persons under bed rest orders</li> <li>○ Persons with acute gastrointestinal complaints or ailments such as gastritis, enteritis, or hepatitis</li> <li>○ Persons under a doctor's care, or those experiencing unusual physical discomfort</li> <li>○ Joint dysfunction due to rheumatism, gout, or other diseases</li> <li>○ Persons with symptoms such as qualm, dizziness, or tinnitus</li> <li>○ Inflammation due to sprain or torn muscles</li> <li>○ Persons with thecitis or suspected thecitis</li> </ul> </li> <li>● <b>Elderly people and persons with weak bones, even if they have no particular disease, should seek medical advice before using the chair.</b></li> <li>● <b>If your condition does not improve after use, consult a doctor or specialist.</b></li> <li>● <b>Do not drop or insert any object into any opening.</b></li> <li>● <b>To disconnect, switch all controls to the off position, then remove plug from outlet.</b></li> </ul>
	<ul style="list-style-type: none"> <li>● <b>Do not use the chair with anything other than an AC 120 V power source.</b> → May result in fire or electric shock.</li> <li>● <b>Do not use with any type of transformer device.</b> → May result in damage to the product or electric shock.</li> <li>● <b>Do not use if the power cord or plug(s) is damaged, or if the wall outlet is loose.</b> → May result in electric shock, short circuit, or combustion.</li> <li>● <b>Be sure to check for tearing on the chair before use. Also, check the other parts for tearing. Stop using the chair immediately if you find even the smallest tear. Unplug the chair and contact the vendor for assistance.</b> → Using the chair with holes or tears in the covering could result in injury or electric shock.</li> <li>● <b>Do not operate under blanket or pillow.</b> → May result in burns, fire, electric shock, or injury to persons.</li> <li>● <b>Do not use outdoors.</b> → May result in burns, fire, electric shock, or injury to persons.</li> <li>● <b>Do not operate where aerosol (spray) products are being used or where oxygen is being administered.</b> → May result in burns, fire, electric sock, or injury to persons.</li> <li>● <b>Keep the cord away from heated surfaces.</b></li> <li>● <b>Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair and the like.</b></li> </ul>

# SAFETY PRECAUTIONS

## Safety Precautions: Before Use

### Caution



Required Action

- **Before plugging the power cord into the wall outlet, make sure that the main POWER switch is set to OFF.**  
→ May result in damage to the product or electric shock.
- **Make sure that the plugs on the power cord are inserted fully into the cord connector and the wall outlet. Do not apply pin, dust or water on the plug.**  
→ May result in electric shock, short circuit, or combustion.
- **Before using, turn on all the switches in order as indicated in this manual to check that the chair is operating properly.**
- **Start off with a lighter massage. Select a massage intensity depends on your physical condition.**  
→ Starting with a deep massage can result in injury. Elderly people and persons with weak bones should be especially careful when selecting massage intensity.
- **Do not sit, stand, or lean on the leg rest unless you are sitting in the chair or getting up from the chair.**  
→ Personal injury or damage to the chair may result if the chair tips over.
- **If you will be using the chair after it has not been used for a long time, first read this Operating Manual again thoroughly and check if the chair operates properly.**



Prohibited

- **Do not use the product in any way not prescribed by the Operating Manual. Do not use the product together with any other therapeutic device or with an electric blanket.**  
→ May result in injury or have a negative effect on health.
- **Do not move the product while someone is in it.**  
→ Accident or injury may result if the product tips over.
- **Do not use with anything other than the accessories supplied with this product.**

## Safety Precautions: During Use

### Warning



Prohibited

- **Do not wrap the power cord around the main body or place the chair on the power cord. Do not scratch, tear, treat, unduly twist, stretch, or bend the power cord.**  
→ The cords may fray or tear, resulting in fire or electric shock.
- **Do not spill water and/or other liquids on the chair or control panel.**  
→ May result in electric shock, short circuit, or combustion.
- **Do not use the chair while the body is wet, or with wet hands.**  
→ Could result in electric shock or damage to the chair.
- **This chair is not intended for use by persons (including children) with reduced physical, sensory or learning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the chair by a person responsible for their safety. Children should be supervised to ensure that they do not play with the chair.**  
→ Could result in accident or injury.

### Caution






Prohibited




- **Do not use the chair unclothed.**
- **Do not use the chair for purposes other than massage, or other than the ordinary uses of a chair. (Do not stand on the main body use as a stepladder etc.)**  
→ Could result in accident or injury.
- **Do not use the chair to massage the head, chest, stomach, and/or joints (knees, elbows).**  
→ May result in nausea or injury.
- **Do not sit wearing hard objects.**  
→ May result in personal injury.
- **Do not operate the product or leave it operating when there is no one in it.**  
→ Could result in accident or injury.
- **Do not use the product with two or more persons at the same time.**  
→ Could result in accident or injury.
- **Do not place anything other than legs and feet on the leg rest.**  
→ May result in personal injury.
- **If your leg or foot comes out of the leg rest during the massage, do not try to force it back in.**  
→ Undue force may result in damage to the mechanism.

# SAFETY PRECAUTIONS






## Safety Precautions: During Use

 <b>Caution</b>	
 Required Action	<ul style="list-style-type: none"> <li>● In an emergency, or if there is any sudden malfunction, push the STOP button on the control panel to make the chair stop immediately.</li> <li>● If you begin to experience marked pain or other unusual physical discomfort while using the product, seek medical advice immediately and discontinue use of the product.</li> <li>● If you feel the massage is too strong, reduce the stimulation according to this Operating Manual. If this solution does not work, stop using and contact the vendor. → May result in injury if you continue use with high stimulation under such circumstances.</li> <li>● Do not fall asleep in the chair during a massage. → May result in personal injury.</li> <li>● Keep massage sessions to no more than 15 minutes in length. → Continual use for extended periods of time can result in over-stimulation, which could cause the opposite of the desired effect, or even injury. Set 30 minutes in total for reference per day and observe appropriate intervals between massage sessions.</li> <li>● Before using, keep person (especially children) or pets off the chair away from the chair. → May result in accident or injury.</li> </ul>
 Unplug the chair	<ul style="list-style-type: none"> <li>● If there is a power outage or the possibility of a power outage, stop the chair immediately and set the main POWER switch to OFF. Unplug the chair from the wall outlet. → Accident or injury may result if the chair is still in operating mode when the power resumes.</li> <li>● If there is a threat of lightning, stop the chair immediately and set the main POWER switch to OFF. Unplug the chair from the wall outlet.</li> <li>● If the chair does not start, or if there is any type of malfunction, to prevent accidents it is best to stop the chair immediately and set the main POWER switch to OFF. Unplug the chair from the wall outlet. Contact the vendor for assistance regarding repairs or maintenance. → Electric shock or fire caused by a short circuit may result if the chair is left plugged in under such circumstances.</li> </ul>

## Safety Precautions: After Use or When Not in Use


 <b>Caution</b>	
 Required Action	<ul style="list-style-type: none"> <li>● Be sure to set the main POWER switch to OFF after every use and to unplug the chair from the wall outlet. → Accident or injury may result if children or pets play with the product.</li> <li>● Wipe the product off after every use, and store it in a place free of excess humidity.</li> <li>● Do not place any object on the product.</li> </ul>
 Unplug the chair	<ul style="list-style-type: none"> <li>● Do not pull on the power cord to unplug the product from the wall outlet; make sure to grasp the plug firmly when unplugging. → May result in electric shock or fire caused by short-circuiting.</li> <li>● When cleaning or caring for the product, make sure to unplug the product from the wall outlet first. Do not plug or unplug the power cord with wet hands.</li> <li>● If you will not be using the product for some time, be sure to unplug it from the wall outlet. → Insulation degradation by dirt or moisture may occur and cause electric shock, or fire caused by short-circuiting may result.</li> </ul>

## Others Precautions


 <b>Warning</b>	
 Prohibited	<ul style="list-style-type: none"> <li>● Do not stand on the chair. Do not sit on any part of the chair other than the seat. → Personal injury or damage to the chair may result if the chair tips over.</li> </ul>
 Do not disassemble	<ul style="list-style-type: none"> <li>● Do not remodel the chair. Anyone other than an authorized service person must not repair or disassemble. Contact the vendor.</li> <li>● For safety reasons, if the power cord is damaged, contact the vendor for replacement.</li> <li>● Do not remodel the chair. → May result in fire, electric shock, or injury.</li> </ul>
 <b>Caution</b>	
 Required Action	<ul style="list-style-type: none"> <li>● If you lend, sell, or otherwise transfer the chair to another party, make sure to include this Operating Manual.</li> </ul>
Note	
<ul style="list-style-type: none"> <li>● Before discarding the product or any of its accessories, consult your municipal government about the proper disposal method.</li> </ul>	

# COMPONENTS & DESCRIPTION

## Chair



## Control Panel



### POWER Lamp

Illuminates to indicate that the power is ON.

### Preset Programs Buttons

Choose from one of four preset programs.

### INTENSITY ADJUSTMENT Button

The air massage intensity can be adjusted to either High or Low.

### STOP Button

The air massage can be cancelled. Press the button to stop the massage immediately in an emergency or if you suspect a malfunction.



Caution

- Do not spill any liquids on the control panel.
- Do not handle the control panel with wet hands.

- Do not use pointed or hard objects to push the buttons on the control panel.  
→ May result in damage to the control panel.
- Do not sit on this panel. Do not put anything on this panel.

# MAIN FEATURES

## Main Features

### From your back to the soles of your feet Air Massage Function

INADA CUBE's uses soothing air to provide shiatsu massage and rotational stretching for your back, hips and seat. Air cells knead your calves, heels, tops of your feet, and the soles of your feet. Shiatsu nodules increase stimulation of your calves and feet, delivering more relaxing, deep therapy through pressure-point stimulation.

### Functional and compact design

Simply fold the back rest of the chair and INADA CUBE instantly becomes a stool or ottoman.

### With state-of-the-art technology Various Course Functions

INADA CUBE offers these four preset programs: FULL, LEGS, STRETCH, and SWING.

## Benefits

As a therapeutic massage device in the general households, the INADA CUBE provides the following benefits:

- Relaxes muscles
- Relives minor muscle aches and pains

# BEFORE USE


- Do not expose to sudden temperature changes. If the room where the chair is to be used is cold, turn the heat up slowly. If the chair is moved from a cold place into a heated room, wait one hour before using the chair.

→ Condensation could occur in the moving parts of the product, which may compromise performance.

## Installation

Keep at least 26 inches behind the chair when installing, so as to be able to install the seat back without obstruction.

- To avoid damaging floor surfaces and to cut down on the noise of the chair, it may be better to install the chair on a carpeted surface or to lay a piece of carpet under the chair.
- Install the chair at least 1 m(3 feet) away from TVs, radios, and/or other audio or video equipment. Interference from these devices can cause problems.



Caution




- Do not install the chair in an area of high humidity, such as a bathroom or sauna.
- Only install on a flat, level surface.
- Do not use in areas of high debris or dust content. Do not use where ambient temperature is 104°F or over.
- Do not install the chair where it can come into contact with direct heat, such as direct sunlight or a radiator.
- Keep away from heat sources, open flame, lit cigarettes, and ashtrays.

# BEFORE USE


## Massage Mode / Stowed Position

### Massage Mode

1. Pull the leg support open until it clicks and locks into place.
2. Pull the seat back slowly while holding the leg support and the seat back as illustrated below.
3. Check to insure that the leg support is stable then carefully sit down.





When fully stowed, the chair can be use as a stool.



### Stowing

1. Use the leg support and seat back to fold into the stowed position is shown below.
2. Press the unlock button at side of the leg support.
3. Replace the leg support to the flat, stowed position.



**Caution** Do not use the chair for purposes other than massage, or other than the ordinary uses of a chair.  
(Do not stand on the main body use as a stepladder etc.)

### Moving

The chair is equipped with wheels to make it easy to move. Turn the Main Power Switch OFF, unplug the power cord from wall outlet and stow the backrest. Hold the leg support as illustrated, and then move the chair on the wheels by lifting, holding up the front of the chair, and gently pushing.




**Caution** Do not move the product with someone is in it.  
 Do not move the product with someone is in it.  
 During moving, watch your toes.  
 During moving, use caution as casters may scratch delicate floors.

**Caution** Be careful not to get hands or objects caught while changing product configuration.  
 Do not place on top of any power cord.

# BEFORE USE

## Checking the Seat Cover

Check for holes or tears at seat back and seat covering.




Warning


- Be sure to check for tearing on the chair before use. Also, check the other parts for tearing. Stop using the chair immediately if you find even the smallest tear. Unplug the chair and contact the warranty provider for assistance.

## Main Power Connection

1. Plug the power cord into the wall outlet.



2. Turn ON the main POWER switch located at the back of the chair.



Warning

- Do not use with any type of transformer device.  
→ May result in damage to the product or electric shock.




Caution

- Before plugging the power cord into the wall outlet, make sure that the main POWER switch is set to OFF.
- Make sure that the plugs on the power cord are inserted fully into the cord connector and the wall outlet.
- Before using the chair for the first time, turn on all the switches in the order indicated in this manual to check that the chair is operating properly.

# INTENSITY ADJUSTMENT

## How to use the Buffer Pad

Your INADA CUBE includes a Buffer Pad that will decrease the intensity of the massage on the soles of your feet. Shiatsu nodules are located under your feet, and the buffer pad can be used to adjust the level of stimulation to suit your personal preference.



Caution





- Be aware of these shiatsu nodules when you sit down or stand up from INADA CUBE.




# OPERATION

## Preset Programs

Four automatic programs provide a variety of different massages to relieve fatigue.

- |   |                |  |
|---|----------------|--|
|  | <b>FULL</b>    | For a comprehensive massage covering your back to the soles of your feet. This session is recommended for relief from general fatigue. |
|  | <b>STRETCH</b> | This session focuses on rotation and rocking motions of hips and low back.   |
|  | <b>LEGS</b>    | This lower body session massages your hips, calves, and feet.  |
|  | <b>SWING</b>   | This session stimulates the hips and buttocks. It's a great overall relaxation massage.  |

## Operation



1. Turn ON the Main Power Switch.
2. Press the desired Preset Program button.  

Press the [Hi/Lo] button to adjust the air intensity. The default setting is "Lo".
3. The program starts, and will run for 15 minutes before stopping.

- While a Preset Program is running, you can push another program button to change the program.

## Stopping in the Middle of a Massage

To stop mid-massage, in an emergency or if there is a malfunction during the operation, press the STOP button. Then turn OFF the Main Power Switch and unplug the power cord from the wall outlet.



- If you feel the massage is too strong, reduce the stimulation according to this Operating Manual.  
→ May result in injury if you continue use with high stimulation under such circumstances.
- Do not sit wearing hard objects.  
→ May result in personal injury.
- If your leg or foot comes out of the leg rest during the massage, do not try to force it back in.
- Keep massage sessions to no more than 15 minutes in length.  
→ Continual use for extended periods of time can result in over-stimulation, which could cause the opposite of the desired effect, or even injury. Set 30 minutes in total for reference per day and keep appropriate intervals between each massage session.

## AFTER USE

Turn OFF the Main Power Switch and unplug the power cord from the wall outlet.



- Be sure to turn OFF the main POWER switch after every use and to unplug the chair from the wall outlet.  
→ Accident or injury may result if children or pets play with the product.
- Wipe the product off after every use, and store it in a place free of excess humidity.
- Do not pull on the power cord to unplug the product from the wall outlet; make sure to grasp the plug firmly when unplugging.
- If you will not be using the product for some time, be sure to unplug it from the wall outlet.

# TROUBLESHOOTING

If there is still a problem after trying the solutions given in this troubleshooting guide, stop using the product immediately so as to prevent accident or injury and unplug the product from the wall outlet. For technical service, please refer to the last page of this booklet.

- Problem** Noise during operation:
- Sound of the air pump
  - Sound of the air massage cushions

**Cause & Solution** These noises are a normal part of the product's operation and do not affect its performance in any way.

- Problem** Control panel buttons do not seem to be activating normally.

**Cause & Solution** Check for the following possible problems:

- ① Check if the power cord is fully plugged into the wall outlet.
- ② Check if the main POWER switch at the rear of the chair is in the ON position.

Anyone other than an authorized serviceperson must not repair or disassemble.

# FAQ

**Q** I spilled water on the product. What should I do?

**A** Immediately turn OFF the main POWER switch located at the rear of the chair, and unplug the chair from the wall outlet. Contact the warranty provider for assistance.

# CARE AND STORAGE

## Care and Maintenance



Caution

- When cleaning or caring for the product, make sure to unplug the product from the wall outlet first. Do not plug or unplug the power cord with wet hands.

### Caring for the Covers

Soak a piece of soft cloth in lukewarm water in which a small amount of mild detergent has been added. Then thoroughly wring it out before cleaning off any spots or grime.



Caution

- The fabric color can change if the product is in contact with vinyl covers or some wall coverings for a prolonged period of time.


### Care for the Footrest Cover and the Buffer Pad

The Footrest Cover and Buffer Pad are machine-washable. Here's how you remove the footrest cover:

1. Pull the footrest cover inward as shown in the picture.



2. Open the zipper and remove the cover.



### Washing instructions...

- Wash separately in case colors bleed.
- Use a fabric softener.
- Air dry in the shade. Do not use a dryer.



Warning

- Never use this chair without the footrest cover. Do not use INADA CUBE if the footrest cover is torn.

### Caring for Other Plastic Parts

Soak a piece of soft cloth in lukewarm water in which a small amount of mild detergent has been added, thoroughly wring it out, and clean off any spots or grime. Then wipe with a dry cloth so that no detergent residue remains.

- Do not use alcohol based or other household cleaning products or abrasives because they can scratch or crack the finish and/or cause color to fade.

### Caring for the Control Panel

Wipe dirt off the control panel with a dry cloth.

- Never use a moistened or wet cloth.  
→ May result in mechanical failure.

## Storage

- After wiping off any dust or dirt, store the product in a place that is relatively dry.
- If you are not going to use the product for some time, put a cover on it that will protect it from dust and dirt.
- Pets can sometimes chew on power cords, so try to keep pets away from the product even when it is being stored.



Caution

- Do not install the product where it can come into contact with direct heat, such as direct sunlight or a radiator.

# CUSTOMER SERVICE AND WARRANTY

## 1. Warranty

Please read the contents carefully and keep it in a safe place. Please complete and mail in the warranty card that was include with your chair.

If the chair is to be used outside the home (in a commercial setting, for example) fees will be charged for repairs.

## 2. Repairs & Service

Before calling for repairs or service, use the Troubleshooting guide on page 17 to see if the problem can be solved easily.

If the problem persists, contact the warranty provider for assistance. Have the following information ready.

- 1) Your name, address, and phone number
  - 2) Product name and model number (see the back cover of this operating manual)
  - 3) Serial No. (see the Serial No. label on the back of the INADA CUBE)
  - 4) Date of purchase
  - 5) Trouble status
  - 6) Proof of purchase
- If the warranty is valid:  
A service charge may apply, subject to the nature of the trouble and/or repair according to the terms of the warranty. See the warranty for full details.
  - If the warranty has expired:  
See the back of this booklet and contact our technical service department.  
If the product can be repaired, we will give you an estimate of the cost to repair the product.
  - Replacement parts stock:  
We will stock service parts (except fabric) of this product for at least six years after terminating the production of this model.  
Service parts refer to those necessary for maintaining the intended performance of the product.
  - For more information on service:  
Please see the back page of this booklet.

# REGULAR MAINTENANCE FOR SAFETY

Regular maintenance will preserve the life of your chair

## Have you experienced any of these problems?

- The chair does not start with the main POWER switch set to ON.
- The power cord has been cracked or scratched.
- The power cord or plug becomes abnormally hot.
- The chair produces a burning smell, abnormal sound or vibration.
- The auto timer function does not work in some cases.
- Other failures.

If you find anything listed above, stop using and unplug the power cord from the outlet to prevent failure or accident. Contact us for maintenance.

## To Use Safely

Fabric of the chair (roller fabric) is consumable, and its service life is 3 years (conditions of use: two 15-minute sessions a day, every day). For safe use, periodic inspection by us is recommended before the end of service life.

- \* Durable years vary depending on conditions of use.
- \* After about five years of use, a professional inspection is recommended even if there are no symptoms or problems.  
For inspection fees, call our technical support number.

# ***SPECIFICATIONS***

---

Name & Model Number	INADA CUBE FML-3000A
Power Supply	AC 120V
Power Consumption	40W
Power Frequency	60Hz
Rated Time	15 minutes
Preset Programs	4 types
Air Pressure	Approx.40kPa
Auto Timer	Approx.15minutes
Dimensions	Approx.19 × Approx.22(39) × Approx.18(28) inch.(W/ D/H) The depth and height figures in () represent the dimension when in the Massage Mode.
Weight	Approx. 49.6lbs.
Outer Covering	Synthetic leather with antibacterial properties Synthetic knit fabric
Accessories	Buffer Pads

Manufactured by FAMILY CO.,LTD.  
2-1-3 Nishi-Miyahara, Yodogawa-ku, Osaka 532-0004, Japan

Made in CHINA  
772294501