

MODERNWATER

# RaPID Assay<sup>®</sup>

## Total BTEX/TPH

RaPID Assay<sup>®</sup> is a rapid field or laboratory enzyme immunoassay method for the analysis of soil and water for remediation, assessment, and industrial testing.

### Test result type

- Quantitative, semi-quantitative or qualitative

### Samples per kit

- 100 test kit (tests up to 80+ samples)  
(Five kit minimum order)

### Assay range

- Soil: 0.9 ppm to 30 ppm as total BTEX standard
- Water: 0.02 to 3.0 ppm as total BTEX standard
- as TPH: range varies based on fuel source

### Sample preparation

- Soil samples require prior extraction using the sample extraction kit (sold separately).
- The sample extraction kit provides materials for 12 soil sample extractions with methanol.

### Sampling time

- Soil extraction time is typically two minutes per sample plus assay run time of approximately 60 minutes.

- Rapid field testing procedure for analysis of soil and water samples
- Quantitative data results with excellent analytical precision
- Results available in approximately 60 minutes
- Training available
- Magnetic particle immunoassay
- Correlated to TPH methods (knowledge of fuel source required)



MODERNWATER



US: +1 302 669 6900  
UK: +44 (0)1483 696 030  
[www.modernwater.com](http://www.modernwater.com)  
[info@modernwater.co.uk](mailto:info@modernwater.co.uk)

## Specificity

The total BTEX/TPH concentrations of the indicated contaminant types in soil samples are expressed below, at each of the three kit calibrator levels, in units comparable to results from GC method 8015 (TPH).

The total BTEX/TPH RaPID assay immunoassay test does not differentiate between BTEX and other petroleum hydrocarbons. The table below shows compounds at the limit of quantitation (LOQ) - an approximate concentration required to yield a positive result at the lowest standard. The IC50 is the concentration required to inhibit one-half of the colour produced by the negative control. It is also used to calculate cross-reactivity values to similar compounds.

TOTAL BTEX/TPH RAPID ASSAY - TOTAL FUEL IN SOIL (PPM)			
Contaminant	Std. 1 equiv.	Std. 2 equiv.	Std. 3 equiv.
Gasoline	10	50	500
Jet A	200	500	2500
JP-4	20	80	750
JP-5	40	100	500
Turbine (Jet) Fuel	10	50	1000
Kerosene	50	150	1250
Artic Diesel	3	10	60
Diesel Fuel	10	40	300
Home Heating Oil	20	70	600
Fuel Oils	10	50	500

BTEX IN SOIL (PPM)		
Contaminant	LOQ	IC50
BTEX*	0.9	6.55
m-Xylene	2.5	18.0
p-Xylene	4.3	31.0
o-Xylene	6.5	47.0
Ethylbenzene	10.8	78.0
Toluene	10.2	74.0
Benzene	70.6	510
Naphthalene	0.82	5.9
1,2,4 - Trimethylbenzene	1.6	11.5
Anthracene	38.8	280
Styrene	3.6	26
Phenanthrene	2.2	16.0
Creosote	6.6	47.6
1,3,5 - Trimethylbenzene	4.8	35.0
Acenaphthene	85.2	620
n-Propylbenzene	6.5	47.0
n-Heptane	99.7	720
n-Hexane	>138	>1000
n-Octane	>138	>1000
n-Nonane	>138	>1000
Cyclohexane	>138	>1000
n-Decane	>138	>1000
Methylene Chloride	>1000	>1000
Trichloroethylene	>1000	>1000

### Test kit components

- Antibody coated magnetic particles for analysis of 100 test tubes
- Zero standard, wash, enzyme conjugate, colour development and stop reagents
- Standards for 0.09, 0.35 and 3.0 ppm as total BTEX
- Kit control as 2.1 ppm as total BTEX
- Disposable test tubes

### Storage & precautions

- Shelf life is typically one year from date of manufacture, with specific kit expiration date information provided on product packaging
- Reagents must be stored at 39° to 46°F (4° to 8°C) when not in use
- Store at ambient temperature 64° to 81°F (18° to 27°C) is acceptable for day of use
- Kits must be brought to 64° to 81°F (18° to 27°C) before use
- Do not expose colour solution to direct sunlight

### Required test materials

- Total BTEX/TPH 100 tube kit
- Sample extraction kit (for soil samples only)
- Total BTEX sample diluent 100ml (as needed to increase range)

### Part #

A00162  
A00165EB  
A00164

### Required test equipment

- RaPID Assay accessory kit which contains RPA-II RaPID analyser
- Magnetic separation rack
- Repeater pipet
- Adjustable volume pipet
- Vortex mixer
- Portable balance
- Digital timer
- Repeater pipet tips
- Adjustable pipet tips

6050100 (p)  
6997010 (r)

6000111  
A00004  
A00008  
A00176  
A00014  
A00131  
A00015  
A00009  
A00013

### Other recommended materials

- Latex gloves
- Liquid and solid waste containers
- Calculator
- Absorbant paper for blotting
- Marking pen

MODERNWATER