


# DOPKE<sup>™</sup> DECON SYSTEM

PureLine

The Dopke<sup>™</sup> Decon System is a bench top sealed chlorine dioxide treatment system for decontaminating PPE like N95 masks, gowns, glasses, shields etc. The system produces and maintains ~1000ppm ClO<sub>2</sub> gas (not a mist, but a "true gas" that obeys ideal gas laws) within the decontamination chamber for just over an hour to provide 6-Log reduction (sterilization) of all surfaces and fibers (includes layered fabrics and filter mask material) inside and out. The system utilizes PureLine's patented PureVista<sup>™</sup>(PV) chlorine dioxide production canisters to safely and easily production gas as needed for each decontamination cycle.

BASE DOPKE™ DECON SYSTEM – Gasketed Transparent Box with PV launching Tray, Sensor Port installed.

The Dopke<sup>™</sup> Decon System's gasketed box is ideal for treating up to 75 N95 masks with chlorine dioxide gas product using PureLine's PureVisa technology. The Decon Box offers a gasketed seal to safeguard operators from gas leaks keeping high levels of ClO<sub>2</sub> gas within the container for optimal decontamination. Durable, tight-fitting latches keep the lid secured to the base, and the clear lid and base allow contents to be easily identified, without opening. The 2 of foot size is ideal for blankets, clothing, masks, and other PPE that need to be decontaminated. The overall assemble dimensions of this item are 22-1/2" L x 16" W x 15-3/4" H (without ClO<sub>2</sub> gas sensor). ClO<sub>2</sub> can be safety released from the system using a vented hood or other active exhaust or a carbon filter recirculation option can be added.



**PureLine** The CIO<sub>2</sub>GÖ<sub>2</sub>

### **OPTIONS**

### LEVEL 1 BASIC

#### DOPKE BOX

A sealed box, where the masks are placed, will hold and decontaminate up to 75 masks at a time. Must be able to ventilate safely. No gas monitor included. Box Dimensions: 22-1/2" L x 16" W x 15-3/4" H.

### LEVEL 2 PREMIUM

#### DOPKE BOX WITH CARBON FILTER DEGASSING UNIT

Contains a recirculation carbon filter system to remove chlorine dioxide from the box prior to opening and removing items. This eliminates the need for ventilation or respirators when the system is opened after a decontamination cycle. Carbon filter is rated for amore than 25 cycles. Replacement cartridges available. Gas sensor not included.

## LEVEL 3 ULTIMATE

DOPKE BOX WITH CIO2 GAS REMOVAL SYSTEM

#### AND CONCENTRATION MONITOR

System with CIO<sub>2</sub> gas sensor and system for data logging USB connection to a PC. PC software included for monitoring and report generation. Capture, record and store all data for up to 100 decontamination cycles. The concentration monitor measures the gas concentration in the Dopke Box for the 75 minute cycle. The sealed carbon filter air circulation system is installed for scrubbing the CIO<sub>2</sub> from the Dopke Box. Prior to opening

the lid CIO<sub>2</sub> levels are reduced to near zero within 15-30 minutes after decontamination is complete. The operator can safely open and remove decontaminated materials without the use of a respirator or exhaust system. Ready to use out of the box.

#### **REQUIRED\*** PV2, EACH WITH LAUNCHER\*

The chlorine dioxide gas dispenser. 5 to a package (5 chlorine dioxide cannisters and 5 launchers). Place 1-2 cannisters per cycle.



#### **OPTIONAL ITEMS:**

#### **DE-GAS NEUTRALIZATION SYSTEM**

Used to neutralize the gas before you open it. Required, if proper ventilation is not available.


#### DATA LOGGING - PC LINK

Ready to Use.

Using a PC you can capture all the time and concentration data for all decontamination cycles.

#### DATA LOGGING - STAND ALONE

Complete stand-alone data logging system. Ethernet based. A PC is not needed for this unit.







PureLine

TM