



Be Right™



MUG Tube

Hach's MUG products (4-methylumbelliferyl-beta-D-glucuronide) provide a safe, prepared, and effective method of confirmatory testing for e.coli bacteria. Using prepared media decreases the risk of contamination and a false positive result. MUG products can be used in the EPA-approved method for E.coli detection, and are available with or without tubes.

Part Number	Description	Quantity	Platform	Material	Platform	USD Price
2182115	Lauryl Tryptose (with MUG) Broth Tubes, concentrated, pk/15	15/pk		Glass		Contact Hach
2217515	Lauryl Tryptose w/ MUG Broth Tubes, Single Strength, pk/15	15/pk		Glass		Contact Hach
2282415	EC MEDIUM/MUG TUBES PK/15	15/pk		Glass		Contact Hach
2471515	EC/MUG without Durham Tubes, Package of 15	15/pk		Glass		Contact Hach