

Accessories

Universal Harness

Attaches to tray to secure low profile plates.

Description	Used on Shakers	Part Number
Universal Harness	Advanced 1000-3 / 1000-15	980192

Universal Harness

**Dimpled Mat**

Designed to hold centrifuge tubes, vials, culture tubes, and micro-tubes securely in place. Mat can easily be removed for cleaning and transporting of tubes from bench to tray.

Description	Used on Shakers	Part Number
Dimpled Mat, 12.75 x 10" (32.4 x 25.4cm)	Advanced 1000RS Rocking Shaker	980535
Dimpled Mat, 11.75 x 8.75" (29.9 x 22.2cm)	Advanced 1000WS Waving Shaker	980542
Dimpled Mat, 10 x 7.5" (25.4 x 19.1cm)	Professional 1000RS Incubating Rocking Shaker	980536
Dimpled Mat, 9.25 x 7.25" (23.5 x 18.4cm)	Professional 1000WS Incubating Waving Shaker	980543
Dimpled Mat, 11.75 x 8.75" (29.9 x 22.2cm)	Advanced 1000-3 / 1000-15	980190
Dimpled Mat, 14 x 11" (35.6 x 27.9cm)	Standard 1000RS Rocking/1000WS Waving Shakers	980TADMAR

Dimpled Mat

**Stacking Tray**

Easily attaches to the units' included tray to add a second tier for higher capacity applications. Second tier tray mounts 3.5" (8.9cm) above lower tray. The tray includes hardware and a rubber mat.

Description	Used on Shakers	Part Number
Stacking Tray, 12.75 x 10" (32.4 x 25.4cm)	Advanced 1000RS Rocking Shaker	980534
Stacking Tray, 11.75 x 8.75" (29.9 x 22.2cm)	Advanced 1000WS Waving Shaker	980541
Stacking Tray, 14 x 11" (35.6 x 27.9cm)	Standard 1000RS Rocking Shaker	980TASTKTRAY

Stacking Tray

**Dilution Cup Tray**

Constructed of type 304 stainless steel. Holds 24 x 28mm dilution vials.

Description	Used on Unit	Part Number
Dilution Rack Tray, 6.9 x 10" (17.5 x 25.4cm)	Advanced 1000-3 / 1000-15	980200

Dilution Cup Tray

**Micro-Tube Rack**

Optional 1.5 to 2mL Micro-Tube Rack attaches to tray to hold up to 70 x 1.5mL or 2mL micro-tubes. Tray can accommodate up to 2 micro-tube racks.

Description	Used on Unit	Part Number
1.5 to 2mL Micro-Tube Rack	Advanced/Professional 1000MP	980191

Micro-Tube Rack

