# **Hood & Clean Bench**



# **Clean Bench**









#### Clean Bench

Advanced type

HEPA Filter

ВС-Н

- Convenient dual controller (outside and inside of chamber).
- Fully compliant with validation requirements.

Clean I	Bench
---------	-------

Basic type

- Pressure gauge at front as standard
- accessory.
  Easy replacement of HEPA filter.

### **PCR Workstation**

- UV blocking about 99.9% with
- transparent acrylic material.
  Essential for sensitive experiment such as PCR.
- **UV Sterilization Cabinet**
- Sterilization using UV lamp. · Count-down timer.

HEPA Filter	HEPA Filter	-
0	0	0
0	0	0
0	0	-
BC-B	PW	UVC

# **Fume Hood**

Filter UV Lamp Fluorescent Lamp Airflow Model











## **Ductless Fume Hood**

- Remove and circulate fume via built-in filter.
- · No need duct. Convenient moving and installation.

## **Extraction Arm Hood**

- Filter box type.
- · Various hood and arms.

## **Molded Fume Hood**

- Table top type with compact design.
- Easy use with connecting duct system.

#### **Molded Observation Fume Hood**

- Easy use with connecting dust system. • Design by transparent windows at the
- front and the rear, making it easy to observe the experiment process from the both side.

Filter	Selectable Activated Carbon, Acid, Formaldehyde, Ammonia/Amines and HEPA Filter.
Model	DLH

Selectable Activated Carbon, Acid, Formaldehyde, Ammonia/Amines and HEPA Filter.

EFS-33 (Dual), EFS-33 (Single)

Selectable Activated Carbon, Acid, Formaldehyde, Ammonia/Amines and HEPA Filter.(EM-33)

PMH-720 + EM-33

Selectable Activated Carbon, Acid, Formaldehyde, Ammonia/Amines and HEPA Filter.(EM-33)

POH-720 + EM-33



Lab Companion's compose Hood and Clean Bench in optimized products protecting experimenters, environments and samples.