

Gloves/Transfer Equipment

Cotton 11714

- Essential tool for proper weight-handling
- Proper protection for preventing hand residue
- Residue can throw off a weight's calibration and certification for accuracy
- Rice Lake recommends using cotton gloves to handle all precision weights

Plastic 11773

- Plastic gloves are an acceptable handling tool to keep dirt and residue off your weights
- Better than handling precision weights bare-handed
- Plastic gloves are economical, disposable, convenient, and help maintain the integrity of your weights



Gloves

Part #	Description	Price
11714	White cloth glove	\$6.00
11773	Vinyl glove	\$2.50



Transfer Equipment

Part #	Description	Price
11711	Forceps, 3 5/8 in 100 g - 1 mg	\$4.00
50035	Forceps, 4 3/4 in	\$37.00
11747	Plastic bar 200 g - 3 kg, screw knob weights	\$17.00
11712	Metal handle 5 kg +, screw knob weights	\$63.00
65640*	Metal forceps	\$17.00
50037	Weight lifter, OIML 2 kg and 5 kg	\$47.00
50000	Tong assembly, 24 cm x 7.5 cm	\$64.00
152414	Tweezers, premium, 140 mm	\$37.00
157109	Tweezers, premium, 220 mm	\$42.00

*Not recommended for use with polished weights.