


Microbiological Products



LuminUltra ATP

LuminUltra ATP microbiological testing solutions provide feedback on contamination sooner for proactive control!

Most Popular Media

EPA-approved m-ColiBlue24 simultaneously detects total coliforms and E. coli in just 24 hours—no confirmation necessary!

Learn More



Whirl-Pak™ Sampling Bags

These bags feature a white, write-on strip for easy sample identification. Mark bags using an ordinary ballpoint pen, pencil, or felt-tip marker. The Dechlor bag meets National Environmental Laboratory Accreditation Conference (NELAC) standard for sampling wastewater treated with 15 mg/L of chlorination agent. Both bags have 100 mL and 4-oz. fill lines (approximate volume).

Prod. No.	Description
2233199	Bag, Sterile, WhirlPak, 207 mL, 100/pk
2075333	Bag, Sterile, Whirl-Pak with Dechlorinating Agent, 177 mL, 100/pk

Sterile Sampling Bottle with Dechlorinating Agent

Hach Modified Colitag Sterile 120 mL Sample Bottle



Prod. No.	Description
8888006	Hach Modified Colitag Sterile Sample Bottles, 120 mL, 100/pk

Sample Vial

Sample vial contains thiosulfate and custody tie. 120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.



Prod. No.	Description
2634001	Sample Vial with Thiosulfate and Custody Tie, 120 mL, 200/pk