

Bath Circulator



Refrigerated & Heating Bath Circulator



Refrigerated & Heating Bath Circulator

Program type

- Programmable temperature control.
- LCD color touch screen.

Working Temp. Range	-35 ~ 150°C -25 ~ 150°C
Temp. Stability at -10°C	± 0.05°C
Bath Volume	5L / 10L / 20L / 30L
Model	RW3-P



Refrigerated & Heating Bath Circulator

General type

- Various temperature range and capacity.
- Smart refrigerating system.

Working Temp. Range	-35 ~ 150°C -25 ~ 150°C
Temp. Stability at -10°C	± 0.05°C
Bath Volume	5L / 10L / 20L / 30L
Model	RW3

Heating Bath Circulator



Heating Bath Circulator

Program type

- Programmable temperature control.
- LCD color touch screen.

Working Temp. Range	Amb. +10-150°C
Temp. Stability at 50°C	± 0.05°C
Bath Volume	5L / 10L / 20L / 30L
Model	CW3-P



Heating Bath Circulator

General type

- Suitable for precise temperature testing above ambient temperature.
- Temperature deviation alarm.

Working Temp. Range	Amb. +10-150°C
Temp. Stability at 50°C	± 0.05°C
Bath Volume	5L / 10L / 20L / 30L
Model	CW3

Programmable Mode

- Setting and saving 10 programs.
- Up to 10 steps can be set per program.
- Set time up to 99 hours and 59 minutes per step.
- Program repetition possible up to 99 times.
- Time set for increasing / decreasing temperature.
- Pre-heating possible before starting program.
- Displays and controls the status of the device through the computer interface.

