

# Hach LBOD Intellical Probe

## A Breakthrough in BOD Measurement!

Hach LBOD Intellical probe uses patented LDO technology to deliver superior confidence in results!

All the same great benefits of the Hach LDO probe, with the built-in features to measure BOD!

- No membranes to replace.
- Calibration in less than one minute.
- Reliable performance and accurate readings time after time.
- Have confidence in your results!

EPA recommends Hach LDO Method 10360 for use in compliance monitoring!\*

*\*Hach's LDO technology has been approved in numerous states for measuring and reporting Dissolved Oxygen (DO) and BOD. Speak with a Hach representative to learn if your state is approved.*



### Order Information

<b>LBOD10101</b>	Intellical LBOD101 Luminescent/Optical Dissolved Oxygen (LDO) Sensor for BOD Measurements, 1 m Cable
<b>2943100</b>	Disposable BOD Bottles, 300 mL, 100/cs
<b>2943900</b>	Plastic BOD Bottle Stoppers, 25/pk
<b>190901</b>	Stopper, ground-glass for BOD bottles
<b>241906</b>	Bottle, BOD, Cap, 6/pk
<b>253335</b>	Nitrification Inhibitor for BOD, Formula 2533™, TCMP, 35 g
<b>253334</b>	Nitrification Inhibitor for BOD, Formula 2533™, TCMP, 500 g
<b>2918700</b>	BOD Seed Inoculum, Polyseed, pk/50 capsules
<b>2514420</b>	ez GGA - No Pipette GGA, Single Dose BOD Standard, pk/20, 2-mL ampules

<b>2484600</b>	Ampule breakers, media ampules for 2-mL PourRite ampules
<b>1416066</b>	BOD Nutrient Buffer Pillows, concentrate to make 300 mL, 50/pk
<b>1486166</b>	BOD Nutrient Buffer Pillows, 3 mL concentrate to make 3L, 50/pk
<b>2436466</b>	BOD Nutrient Buffer Pillows, 4 mL concentrate to make 4 L, 50/pk
<b>1486266</b>	BOD Nutrient Buffer Pillows, 6 mL concentrate to make 6 L, 50/pk
<b>1486398</b>	BOD Nutrient Buffer Pillows, concentrate to make 19 L, 25/pk

Learn More

